

Esophageal Temperature Monitoring System

# ***EsoMeter***

With wide coverage


### EsoMeter Catheter

Pentapolar / 7F

| Model No. | Electrode Spacing (mm) | Usable Length (mm) |
|-----------|------------------------|--------------------|
| EC5547    | 5-5-5...               | 800                |
| EC5567    | 5-5-5...               | 810                |

Catheter Material

| Electrodes            | Shaft |
|-----------------------|-------|
| Stainless Steel (SUS) | PEBA  |



### EsoMeter Monitor

| Model No. | Monitoring Temperature Range (°C) | Warning Temperature Settings (In Steps of 0.5°C) |                              | Size (mm)<br>[Width x Height x Depth] | Weight (kg) | Power Supply   |
|-----------|-----------------------------------|--|------------------------------|---------------------------------------|-------------|----------------|
|           |                                   | Upper Temperature Limit (°C)                     | Lower Temperature Limit (°C) |                                       |             |                |
| EM-1      | 5~75                              | 36~50  | 5~36                         | 180×160×80                            | 1.7         | AC100V 50/60Hz |

Accessories

| Model No. | Product Name                 |
|-----------|------------------------------|
| EM502     | ECG Cable                    |
| EM507     | Power Cable                  |
| EM510     | Connection Cable (2.8m)      |
| EM512     | Potential Equalization Cable |

| Model No. | Product Name      |
|-----------|-------------------|
| CP509*    | IV Pole connector |
| EM511*    | USB Cable         |
| EM513*    | Check Connector   |

\* Not included as a medical device

Additional Accessories

| Model No. | Product Name            |
|-----------|-------------------------|
| EM510-4** | Connection Cable (4.0m) |
| EM510-5** | Connection Cable (5.0m) |

\*\* Custom-made product: Please contact our local sales representative for further information.

Some models may not be approved in your country. Please contact our local representatives for further information.

Japan Lifeline Co., Ltd.

2-2-20, Higashishinagawa, Shinagawa-ku, Tokyo 140-0002 Japan  
<https://www.japanlifeline.com>